



PETSTAT[®]

Advanced Intelligence for medical imaging

PET/CT Medical imaging viewer

PETSTAT is a PET/CT medical imaging viewer, which focuses on cancer diagnosis, clinical trials and research by PET imaging. PETSTAT has intuitive and easy-to-use GUI which integrates various functionalities such as image analysis, volumetric analysis, histogram analysis, 3D graphics.



●PERCIST support


>Displays SUV bw, SUV bsa, SUV lbm (SUL),Peak

●RECIST support

>Lesion manager lists up measures, ROIs, marks and sums up those parameters

●PET/CT fusion image


>PET/CT fusion image enables quick and intuitive diagnosis.



PET/CT fusion image and Sphere ROI analysis

●Sphere ROI analysis

>Displays SUV max, mean, TLG (Total Lesion Glycolysis) , Peak within a sphere ROI.
>Displays SUV histogram and SUV-cumulative volume and SUV-TLG graph and outputs by CSV format.



●MTV manager

>Semi-automatic MTV (Metabolic Tumor Volume) creation and listing up MTVs with parameters such as TLG , volume, SUV Max, SUV Peak, T/N ratio, Maximum diameter , and summing up these quantities.

>Displays MTVs on MPI/VR views.


>Outputs MTVs by CSV format.

Name	Co...	TLG[g]	Volume...	Thr...	Max	Peak	Me...	T/N r...	LD[mm]	SD[mm]	LDxSD[mm...	Position	Est part
MTV487		3911.15	1999.20	1.48	7.29	4.68	2.80	6.28	171.64	136.77	29474.58	126.15	Brain
MTV381		437.30	191.84	1.48	5.11	4.17	2.28	4.40	98.81	68.00	6718.94	-204.55	
MTV048		337.49	130.13	1.48	4.29	3.57	2.59	3.70	75.18	63.46	4771.23	-581.97	
MTV078		100.67	50.55	1.48	3.85	3.18	1.99	3.32	97.05	37.62	3651.06	-372.81	
MTV080		74.80	38.04	1.48	4.21	3.12	1.97	3.62	84.06	37.60	3160.76	-369.09	
MTV491		36.65	11.31	1.48	10.91	6.27	3.24	9.39	48.71	15.54	756.87	167.06	
MTV091		24.94	12.99	1.48	3.31	2.62	1.92	2.85	43.06	28.76	1289.53	-307.16	
MTV455		20.58	11.77	1.48	2.37	2.03	1.75	2.04	79.46	19.06	1514.34	-90.60	
MTV456		17.09	10.53	1.48	2.27	1.79	1.70	1.95	68.79	21.87	1504.68	48.19	
MTV053		14.85	7.89	1.48	2.94	2.07	1.88	2.53	64.70	15.42	987.32	-483.94	
MTV074		13.25	8.08	1.48	2.18	1.68	1.64	1.88	64.00	18.22	1165.98	-420.24	

●3D graphics

>MIP(Maximum Intensity Projection) which is a prerequisite for PET imaging.

>Volume Rendering which is helpful in intuitive understanding of tumor locations.



●PACS

>Supports DICOM network communication.

●Data output

>Histogram or MTVs can be output to CSV / Text files.

>Stores parameters of viewer sessions for statistical analyses of new indices or effect of cancer therapies.

●Windows platform

>Available on laptop PCs.

>More comfortably available on high-end PC workstations.

●Meeting the needs of cancer frontier

>Reflects the options of the central hospital of National Cancer Center(Japan).

Sales and Development



Ad-In Research, Inc.
<http://www.adin.co.jp>

Address: 2F, Kioi Park building, Kioi-chyo 3-6, Chiyoda-ku, Tokyo, 102-0094 , Japan
Phone: +81-3-3288-7835 Fax: +81-3-3288-7334
E-mail: petstat@adin.co.jp URL: http://www.adin.co.jp/service/detail_9.html